

## Tutorial

1. Download the code from Github. You will get an Intelligent-quantifying-RGCs.zip file which you need to unzip.
2. Download the files in cloud disk that need to be unzipped. The name of the downloaded file is yolov5\_cpu.


Google Link:

<https://drive.google.com/file/d/1yOEsBvil6KEdZFa5ENQxB667uKyEdnKq/view?usp=sharing>


Baidu Link:

<https://pan.baidu.com/s/1lccg1OVbeudsp2VtnqxWZg>  
Extraction code : g44k

3. Find the “Config.ini” file in step 1 and change the file content to your yolov5\_cpu path.
4. By clicking “userinterface.exe” opens the graphical interface. The loading speed is slow when it is opened for the first time. (please run as administrator)



5. After opening the software, you need to check whether the pmse\_plus.pt model file exists in the root directory of drive C. If not, you need to copy it. (The pmse\_plus.pt model file is located in the weights folder of the yolov5 folder, which can be found in the extracted Intelligent-quantifying-RGCs.)



6. Prepare the pictures to be detected. The software only allows one or five pictures to be read at a time. Otherwise, an error will be reported. And it should be noted that the five pictures were taken at different Z-positions in the same field of view.

